



THE CHOICE IS YOURS



BEGO SEMADOS® SC/SCX IMPLANTS
QUICK & FUNCTIONAL · UNIVERSAL & SUCCESSFUL
RELIABLE & PROVEN · VARIED & SOPHISTICATED
100% GERMAN QUALITY PRODUCT
SIMPLE & INTUITIVE · QUICK & FUNCTIONAL
UNIVERSAL & SUCCESSFUL · RELIABLE & PROVEN
VARIED & SOPHISTICATED · SIMPLE & INTUITIVE

Your tried and true system – familiar yet new

Tailored to your patient down to the finest detail. Various implant systems, supplemented by individual CAD/CAM fabricated prosthetic components and surgical guides as well as the BEGO Multi^{Plus} concept – a system that provides patients with a perfectly fitting prosthesis in just one day. The exceptionally high technical standards ensure the highest level of satisfaction – for both dental implantologists and patients.

Your benefits

- SC/SCX implants – preservation of the crestal bone thanks to the integrated Platform Switch
- Efficient – just one surgical tray needed
- Well thought out – cover screw included
- Freedom of choice based on preference and indication – machined implant shoulder (SC implants) or micro-structured shoulder (SCX implants)
- Optional – guided/template-guided preparation of BEGO Semados® SC/ SCX implants Ø 3.25–4.5 using the BEGO Guide System
- For the most stringent quality requirements – 100% German development and production

Technical details

- Gap-free fit of the abutments thanks to the internal taper
- Self-tapping thread
- Manufactured from cold-worked grade 4 pure titanium with ultrapure homogeneous TiPure^{Plus} surface
- Diameters:
3.25 · 3.75 · 4.1 · 4.5 · 5.5 mm
- Lengths:
7 · 8.5 · 10 · 11.5 · 13 · 15 mm
- Diameter 5.5 mm in length 7 mm if there is not much vertical bone height available

PRODUCT OVERVIEW

Distributor contact data

BEGO Semados® SC Implants

Art. no	Colour code	Product
BEGO Semados® SC Implants 3,25		
58160	■	SC 3.25 L8.5
58161	■	SC 3.25 L10
58162	■	SC 3.25 L11.5
58163	■	SC 3.25 L13
58164	■	SC 3.25 L15
BEGO Semados® SC Implants 3,75		
58165	■	SC 3.75 L8.5
58166	■	SC 3.75 L10
58167	■	SC 3.75 L11.5
58168	■	SC 3.75 L13
58169	■	SC 3.75 L15
BEGO Semados® SC Implants 4,1		
58170	■	SC 4.1 L7
58171	■	SC 4.1 L8.5
58172	■	SC 4.1 L10
58173	■	SC 4.1 L11.5
58174	■	SC 4.1 L13
58175	■	SC 4.1 L15
BEGO Semados® SC Implants 4,5		
58176	■	SC 4.5 L7
58177	■	SC 4.5 L8.5
58178	■	SC 4.5 L10
58179	■	SC 4.5 L11.5
58180	■	SC 4.5 L13
58181	■	SC 4.5 L15
BEGO Semados® SC Implants 5,5		
58182	■	SC 5.5 L7
58183	■	SC 5.5 L8.5
58184	■	SC 5.5 L10
58185	■	SC 5.5 L11.5
58186	■	SC 5.5 L13

BEGO Semados® SCX Implants

Art. no	Colour code	Product
BEGO Semados® SCX Implants 3,25		
58260	■	SCX 3.25 L8.5
58261	■	SCX 3.25 L10
58262	■	SCX 3.25 L11.5
58263	■	SCX 3.25 L13
58264	■	SCX 3.25 L15
BEGO Semados® SCX Implants 3,75		
58265	■	SCX 3.75 L8.5
58266	■	SCX 3.75 L10
58267	■	SCX 3.75 L11.5
58268	■	SCX 3.75 L13
58269	■	SCX 3.75 L15
BEGO Semados® SCX Implants 4,1		
58270	■	SCX 4.1 L7
58271	■	SCX 4.1 L8.5
58272	■	SCX 4.1 L10
58273	■	SCX 4.1 L11.5
58274	■	SCX 4.1 L13
58275	■	SCX 4.1 L15
BEGO Semados® SCX Implants 4,5		
58276	■	SCX 4.5 L7
58277	■	SCX 4.5 L8.5
58278	■	SCX 4.5 L10
58279	■	SCX 4.5 L11.5
58280	■	SCX 4.5 L13
58281	■	SCX 4.5 L15
BEGO Semados® SCX Implants 5,5		
58282	■	SCX 5.5 L7
58283	■	SCX 5.5 L8.5
58284	■	SCX 5.5 L10
58285	■	SCX 5.5 L11.5
58286	■	SCX 5.5 L13